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**Description:**

Endotracheal tube for oral and nasal intubation, for patient respiration.  
Manufactured from pure polyvinyl chloride (PVC).

Each individual tube has the following product properties:

- Thermoreactive material, the firmness of which adapts individually to the respiratory tract of the patient in order to avoid or minimise traumatisation.
- Radiologically positive, blue contrast stripe in the lumen of the tube enabling radiological control of the tube's position.
- All tube packs contain a 15 mm connector corresponding to ISO 5356
- All blockable endotracheal tubes are equipped with a stable low-pressure cuff and have an external pilot balloon with one-way Luer lock valve.
- All tubes are equipped with a Murphy eye.
- For depth control a centimetre scale is printed on the tube stating the distance from the distal end of the tube.
- To provide shape, the tubes contain a guide wire to simplify intubation.

**Indications:**

Endotracheal tubes are suitable for oral and nasal intubation of patients for controlled or assisted respiration.

**Contraindications:**

Use of an endotracheal tube is contraindicated if interventions with active or passive laser radiation or electrical radiofrequency currents are planned or can be expected in the direct spatial vicinity of the tube.

Contact with laser radiation or RF surgery can lead to dangerous complications or extreme thermal injury, especially when oxygen is used. Toxins could also be released in melting reactions of the material!

**Warnings:**

1. When shortening the tube with scissors, make sure that the supply line for the cuff running parallel to the tube lumen is not destroyed. Otherwise the cuff could deflate and leakages would occur.
2. Use no viscous lubricants for the lumen of the tube!
3. Deflate the cuff completely before extubation or correcting the position to prevent traumatisation of the trachea.
4. If the tube should be made slippery before the intubation, make sure that the lumen of the tube is not constricted.
5. When anaesthetising gels or lubricants are used, both partial and complete obstruction of the tube have been observed.
6. Caution: Do not overinflate the cuff of the tube! Inflating the cuff too strongly can lead to damage to the tracheal wall, to slipping of the tube and to obstruction of the respiratory tract.
7. Do not estimate the filling volume of the cuff according to the volume of air with which the cuff was filled. Also not by the resistance that you feel in the syringe during the cuff inflating process.
8. During anaesthesia gases can diffuse into the cuff and lead to a rise in pressure inside the cuff!

### Precautionary measures:

1. Check the safety of the entire system and all its connections regularly for firm hold and correct location.
2. The respiration air or the respiration gases should be sufficiently moistened to reduce encrustations of the lumen of the tube and to prevent mucosal lesions!
3. Safe use of the tube should be guaranteed by regular aspiration. Routine checks of the lumen and position of the tube guarantee functional use over the required period!
4. Change or remove the tube at the latest every 30 days!
5. Check the pressure inside the cuff at regular intervals and adapt it if necessary in the case of deviations.
6. The syringe for inflating the cuff must be clean and free of foreign bodies. Remove it from the valve immediately after inflating the cuff!
7. Avoid contact between the tube and sharp objects so that the cuff is not destroyed. If the cuff is defective, do not use the tube or re-intubate the patient and dispose of the defective tube!
8. The supply line to the cuff could hinder MRT scans. Make sure that the valve does not lie in the area to be scanned.
9. Changing the position of the tube with the cuff inflated is not allowed since this can cause injuries!
10. Excessive movement or positioning of head or neck with the tube in place should occur only with applied tube strengthening or when a anti-kink protection is used, so that there is no twisting and consequently obstruction of the tube.
11. The tube should be secured and sufficiently fixed to prevent dislocation.
12. Hygienic standards must be observed in the use of endotracheal tubes.
13. After the patient is repositioned, check the correct position of the tube!
14. Use a suitable bite block.

### Hints for use:

Intubation and extubation must be according to the applicable medical regulations. All applicable instructions with regard to size, shape and type of the tube must be observed so that the correspondingly suitable product is selected!

1. Check the function and tightness of the tube as well as of the cuff before use.
2. Absolutely fit the 15 mm connector before the intubation.
3. Intubate the patient nasally or orally according to the applicable regulations.
4. The depth mark on the outside of the tube is an aid.
5. When positioned correctly, the guide wire must be removed!
6. Inflate the cuff minimally until respiration free of leakage is possible.
7. Check the position of the tube by auscultation or bronchoscopy.
8. Deflate the cuff completely before extubation.
9. Do not use if the packaging is damaged or opened.
10. Protect the product against excessive temperatures and UV radiation.
11. For use one time only! Do not re-sterilise! Do not use if the expiry date is exceeded!
12. The device is for single use only. Reuse risks spreading infectious bacteria. Sterilization and washing (for reuse) may impair performance of the device and causes harm for the patient.
13. The device contains DEHP and is a specific potential risk for the treatment of pregnant or nursing women and children. Therefore the device should only be used for the treatment of these patients after a careful and individual risk-benefit analysis.

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P. J. Dahlhausen & Co. GmbH