---

**ultraMEDIC<sup>®</sup>**

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---

# Contents

	<i>Page</i>
<b>1. Product description</b> .....	<b>2</b>
<b>2. Important information</b> .....	<b>10</b>
2.1 Liability and warranty .....	10
2.2 User group .....	10
<b>3. Information on the safety and protection of the product</b> .....	<b>10</b>
<b>4. Intended use</b> .....	<b>11</b>
4.1 Cautions .....	11
4.2 Use .....	11
<b>5. Cleaning and care</b> .....	<b>12</b>
<b>6. Maintenance</b> .....	<b>12</b>
<b>7. Storage</b> .....	<b>13</b>
<b>8. Service life</b> .....	<b>13</b>
<b>9. Optional accessories</b> .....	<b>13</b>

---

## 2. Important information

Please read the operating manual and safety instructions carefully before you start using the S-BOON. Keep the operating manual in a safe place to avoid operating errors.

Should you require further information or encounter difficulties in its use or maintenance, please contact your local dealer or ultraMEDIC GmbH directly. Always quote the LOT-number to identify the S-BOON.

### 2.1 Liability and warranty

The manufacturer grants a warranty period of 24 months on this product. The warranty only covers the costs of the spare parts. Any incurred shipping or wage costs are not subject to warranty.

Please check the shipment for damage and completeness immediately upon opening. If damage is found, this must be reported within two working days, packaging damage must be reported immediately. After this period, no guarantee or warranty claims for transport damage can be asserted. No responsibility is taken for improper assembly.

### 2.2 User group

This product may only be used by competent, trained and qualified personnel, and the principles of use must be observed. A prerequisite for use is theoretical and practical experience in handling the spineboard S-BOON.

## 3. Information on the safety and protection of the product

1. Before and after every use, check that the components and safety belts are intact and complete.
2. Make sure that the safety belts and components are assembled correctly.
3. Before using the S-BOON for the first time, conduct a test run to familiarize yourself with the product.
4. In the event of damage or suspected damage to the safety belts or

---

S-BOON remove the product from operation immediately until any occurring irregularities have been remedied.

## 4. Intended use

The S-BOON is a medical product for the immobilization of patients with severe head and backbone injuries when they have to be lifted, repositioned or transported.

The "S-BOON" spineboard can also be used for X-ray, MRI and CT scans.

### 4.1 Cautions

1. A specialist medical assessment must be conducted before the injured party is moved, lifted or transported.
2. To ensure the patient's safety, at least 2 users are required to operate the S-BOON.
3. The injured party must be immobilised on the S-BOON during every use. Insufficient or no restraint may lead to severe patient or rescuer injuries.
4. Serious incidents that have occurred in connection with the product must be reported by the user and / or patient to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

### 4.2 Use

1. The S-BOON is applied depending on the injury and local conditions.
2. If the injured person is shorter than 1.20 m, we recommend using the straps SAN-0281 or SAN-0286.
3. After the patient is placed on the spineboard, their head can be fixed with the ultraHEADLOCK I or II.
4. Fixing with the straps is done
  1. at the height of the chest, arms can also be fixed
  2. at the height of the pelvis and
  3. below the patella.

---

Fasten all 3 patient safety belts. The locking of the buckles must always audibly click into place. To be on the safe side, a quick check should be carried out by pulling the appropriate point. Make sure that the patient's breathing is not hindered by pulling the belts and that the belts do not put pressure on the possible injuries.

If the patient is not yet lying securely on the S-BOON, additional straps can be used.

**5.** If the fixing on the spineboard is done using the strap spider system (ultraSTRAPS), the Y-shaped part of the system is positioned towards the shoulder, the rest of the strap is placed over the patient along the middle. The patient is strapped loosely from head to foot. Thereafter, the Velcro straps are tightened in the same sequence till the patient is securely bound on the spineboard.

**6.** Check straps and fasteners before lifting and transporting the patient. Ideally, the transport takes place with 4 to 6 users, who can each reach into the recesses provided (recessed grips).

## **5. Cleaning and care**

Use mild soap and warm water (40°C) to clean the S-BOON.

For disinfection, you can use commercially available, alcohol-based agents. Other solutions or cleaning chemicals may damage the material of the straps, carabiner hooks and clasps.

Allow the wet transportation suspension to dry naturally and keep it away from direct heat exposure.

## **6. Maintenance**

In order to ensure the safety of the injured party and the user, we recommend that you have the S-BOON visually inspected by an expert at least once a year. Appointed auditors must be demonstrably adequately qualified.

For users in Germany, we recommend consideration of the following regulations and guidelines:

DGUV Grundsatz 305-002

DGUV Regel 112-198/199

MPBetreibV

---

**The following parts must be checked and the following defects must be excluded:**

**Base plate:** deformation, tears, breaks, loosened screws, broken parts

**Closures:** impaired function, deformation, torsion, tears and breaks

**Straps:** Cutting and abrasion points, damaged seams, tissue-damaging stains

**Metal parts:** Rust and oxidation

**other damages**

In the event of damage or suspected damage, remove the product from operation immediately.

## 7. Storage

Store the S-BOON in a dry and dark place. Avoid direct sunlight and direct heat radiation from heat sources. These factors could damage the material.

## 8. Service life

The discard criteria depends on the product, its frequency of use and the external conditions of use.

The service life of the S-BOON depends on the individual conditions of use. Depending on the frequency and intensity of use of the transportation suspension, certain wear and tear is normal. Appropriate use is assumed here.

The service life of the optional accessories can be found in the associated operating instructions.

## 9. Optional accessories

SAN-0284	ultraHEAD-LOCK I	Head restraining set
SAN-0285	ultraHEAD-LOCK II	Head restraining set
SAN-0286	ultraSTRAPS Xtreme	Safety strap system with clip locks
SAN-0281	ultraSTRAPS	Safety strap system with Velcro fastener
SAN-0280-T	ultraSPINE BAG	Storage bag
SAN-0084-S	ultraBELT	Replacement safety belt
SAN-0084-K	ultraBELT Set	3 Safety belts + bag

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